



don whitley
scientific

excellence in microbiology

Whitley Workstations Whitley Workstation Accessories



Microscope & Bracket Assembly



Anaerobic Indicator Pump

To complement the Whitley Workstation range of anaerobic and microaerobic workstations, we have developed a range of accessories to enhance your working practices



20 Plate Adaptor Kit



Basket Tag Set

Microscope & Bracket Assembly

A02580

A long focal length microscope will enable the detailed inspection of colonies during early growth periods. Supplied complete with all fixtures and fittings.

Basket Tag Set

A02750

Enables users to colour code individual stacks of Petri dishes. A coloured tag can be attached to any of the Whitley Petri dish wire racks. The re-usable tags are manufactured from polypropylene for durability. Notes can be written using either water soluble or permanent markers.

A set comprises:

- 20 x yellow
- 20 x blue
- 20 x red
- 20 x green
- 20 x white
- 20 x orange
- 20 x grey

20 Plate Adaptor Kit

A02780

Increases the capacity of the standard Whitley Workstation porthole from 10 to 20 Petri dishes. The kit comprises:

- 1 x right porthole attachment, 1 x left porthole attachment
- 2 x thumb nuts
- 2 x 20 plate Petri dish holders

Anaerobic Indicator Pump

A03002 (for use only with the DG250 workstation)

A05019 (for use with any 500/1000 Whitley workstation)

Bubbles internal atmosphere through anaerobic indicator solution to provide an indication of anaerobiosis.

In the interests of a policy of continuous product improvement the company reserves the right to alter product specifications without prior notice. All rights reserved.
© 2006 Don Whitley Scientific Limited.

Don Whitley Scientific Limited

14 Otley Road, Shipley, West Yorkshire, BD17 7SE, England.

Telephone: +44 (0)1274 595728 Fax: +44 (0)1274 531197

Website: www.dwscientific.co.uk Email: info@dwscientific.co.uk